

Outer Leather

Nappa standard L R

Nappa impregnated L R

Buckskin L R

Deerskin L R

Lining

Silk L R

Wool L R

Lambskin, thin L R

Lambskin, thick L R

Deer skin L R

Unlined L R

Colour

Black

Bordeaux

Camel

Brown


gray

Types of fastener


 L R

 L R

 L R



 L R

 L R

 L R

Facing on the palm L R

Quantity



1. Length of all fingers
2. Glove length until the wrist.
3. Fist closure perimeter
4. Glove perimeter at the wrist
5. Glove perimeter at the end of the glove
6. Length between wrist and end of the glove
7. Perimeter

	L	R
thumb	-----	-----
forefinger	-----	-----
middle finger	-----	-----
ring finger	-----	-----
little finger	-----	-----

Please always attach this data sheet to the carefully drawn original. Please also always draw and measure the healthy hand.

Order Estimate for the costs


Company

Contact


Commission

Profession


Taking measurements correctly = optimal fit



1. If possible, always measure both hands
2. Please sketch both hands (provides us with the original shape of the disabled hand)
3. If possible, always sketch the hand straight with the middle finger forming a straight line with the forearm
4. Do not spread the fingers too widely apart



1. If the finger(s) are severely bent,
2. please roll up a thin piece of cardboard,
3. place it in the palm,
4. draw the outline on it and
5. measure the finger length.



1. If there is a protrusion on the back of the hand,
2. please position the hand upright on its side,
3. draw around the hand on its side and
4. measure the protrusion and plot its location.
5. Please still draw the outline of the hand laid flat on the paper.

- For major deformities, please attach pictures.
- If possible, fabrication of a plaster model in neutral position.
- Outline drawing

TIP: Gloves that already fit are highly suitable as patterns, please enclose them.

Please fill out both sides of the datasheet!